



	★ GOOD ★	★★ BETTER ★★				★★★ BEST ★★★				BARIATRIC	
Item #	TR37E-SV, TR39E-SV	ATC17-BK, ATC17-BL, ATC17-RD, ATC19-BK, ATC19-BL, ATC19-RD	RTLFW19RW-RD	EXP19LTBL, EXP19LTRD	TS19	TC005GY	TR18, TR20	SL18	DFL19-BL, DFL19-BLK, DFL19-RD	BTR20-B, BTR20-R, BTR22-B, BTR22-R	ATC22-R
Product Name	Steel Transport Chair	Aluminum Transport Chair	Deluxe Fly-Weight Aluminum Transport Chair	Lightweight Expedition Aluminum Transport Chair	TranSport Aluminum Transport Chair	Travelite Transport Chair w/ Bag	Poly-Fly High Strength, Lightweight Wheelchair/Flyweight Transport Chair Combo	Super Light, Folding Transport Chair w/ Bag	Fly-Lite Aluminum Transport Chair	Bariatric Steel Transport Chair	22" Bariatric Aluminum Transport Chair
Sales Channel											
Key Differences	A steel transport chair available in a 17" and 19" wide seat option.	An aluminum transport chair available in a 17" and 19" wide seat option with multiple colors available.	This lightweight transport chair that folds and compacts for easy storage packaged in a retail box.	Lightweight and a 300lb weight limit. This chair has a comfortable seat width.	This dual color blue transport chair has a sporty look, is easily foldable and packaged in a retail box.	With the lightweight design, this transport chair folds small enough to fit inside it's own carry bag with wheels for easy mobility. It is conveniently packaged in a retail box.	When it comes to unique and different, Poly-Fly is hard to beat. It can be effortlessly transformed from Wheelchair to Transport Chair. The arms are easily removable for easy patient transfer.	Simple folding flip-back design makes using this super light transport chair a dream. It comes with it's own carry bag for storage and transport.	Lightweight and stylish design are the apex of the Fly-Lite transport chair line.	Steel Bariatric Transport Chair with cross brace design folds easily and is offered in 2 seat widths and 2 colors.	Aluminum Bariatric Transport Chair with cross brace design folds easily and is lightweight.
Overall Dimensions Open (W x H)	TR37: 21.5" x 39" TR39: 23.5" x 39"	ATC17: 20" x 39.5" ATC19: 22" x 39.5"	22" x 36"	22" x 38"	24.5" x 36"	22.5" x 37.5"	TR18: 23.5" x 36" 26.5" x 36" (with mag-wheels) TR20: 26.25" x 36" 29" x 36" (with mag-wheels)	23" x 37"	22" x 36"	BTR20: 24.5" x 37" BTR22: 26.5" x 37"	30" x 37.5"
Overall Dimensions Closed (W x H)	TR37: 8" x 39" TR39: 10.75" x 39"	ATC17: 9" x 39.5" ATC19: 9" x 39.5"	9" x 36"	10.75" x 38"	10" x 36"	11.5" x 37.5"	TR18: 12" x 36" TR20: 14" x 36"	9" x 37"	9" x 36"	BTR20: 9" x 37" BTR22: 10" x 37"	12" x 37.5"
Seat-To-Floor Height	19"	19"	19"	19"	18"	19"	19.5"	20"	18"	19.5"	21"
Back Height	18"	18"	18"	18"	16"	17.5"	16"	17"	18"	16"	16"
Seat Dimensions (W x D)	TR37: 17" x 16" TR39: 19" x 16"	ATC17: 17" x 16" ATC19: 19" x 16"	17" x 16"	19" x 16"	19" x 16"	18" x 18"	TR18: 18" x 16" TR20: 20" x 16"	18" x 17"	19" x 15.5"	BTR20: 20" x 18" BTR22: 22" x 18"	22" x 18"
Overall Length (with riggings)	33"	33"	33"	38"	38"	33"	42"	30"	33"	38"	41"
Weight	26 lbs	19 lbs	19 lbs	19 lbs	17.5 lbs (Without Footrests) 20 lbs (With Footrests)	20.5 lbs	34 lbs 36 lbs	19 lbs	18.9 lbs	45 lbs 49 lbs	33 lbs
Carton Shipping Weight	31 lbs	24 lbs	24 lbs	24 lbs	19.5 lbs	23.2 lbs	39 lbs 41 lbs	24 lbs	23.8 lbs	50 lbs 54 lbs	39 lbs
Weight Limit	250 lbs	300 lbs	300 lbs	300 lbs	250 lbs	250 lbs	250 lbs	250 lbs	300 lb.	450 lbs	450 lbs
Retail Box	--	--	✓	--	✓	✓	--	--	--	--	--
Color	Silver Vein Frame, Black Upholstery	Red, Blue or Black Frame and Black Upholstery	Red Frame and Black Upholstery	Red or Blue Frame and Black Upholstery	Blue Upholstery, Blue Frame	Silver Frame and Black Upholstery	Silver Frame and Black Upholstery	Silver Frame and Blue Upholstery	Black or Red Frame and Coordinated Upholstery	Red or Blue Frame and Black Upholstery	Red Frame and Black Upholstery
HPCS	E1038	E1038	E1038	E1038	E1038	--	E1038 or K0004+E0951+E0971+E0978	E1038	E1038	E1039	E1039
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime

LEGEND



E-COMMERCE
Available to order online



RETAIL
Available in large retail chain stores



PROVIDER
Available in independent DME provider locations

GOOD ★

Good models are the base within the class, comparable to the competition of the same level. These meet the basic requirements with less refinements and minimal extras vs. the "BETTER" and "BEST" classifications. They tend to be entry level in the market with good value.

BETTER ★★

Mid-level models within the class, offer more features and benefits to give the consumer a better user experience. These tend to be lighter weight, have more color options available and may have a higher user weight limit.

BEST ★★★

Best models combine similar design features from our "Good" and "Better" models but offer additional features that pushes them to the top ranking. They may fold very small, may be very light or could come with a separate carry bag. They may have unique features, such as the Poly-Fly that is both a Transport Chair and Wheelchair. They tend to include the desirable features for the customer and work to give them the experience they are looking for.