



## SYNERGY AIR ELITE THERAPY SURFACE

The **Synergy** Air Elite Surface is a constant low pressure redistribution surface with microclimate management (LAL)\* and optional alternating pressure designed to help prevent and treat pressure injuries.

Based on a published JWOCN Support Surface Algorithm, a reactive Constant Low Pressure (CLP) support surface with microclimate management (LAL)\* is suggested for patients with moderate to high risk of pressure injuries, limited mobility, excessive moisture and persons with multiple Stage II-IV or sufficiently large Stage III-IV pressure ulcers on the trunk or pelvis involving more than 1 available turning surface.<sup>1</sup>



\*(LAL) Low Air Loss



### FEATURES

- MICROCLIMATE MANAGEMENT  
Continuous airflow removes excess heat and moisture to help keep patient's skin cool and dry.
- ALTERNATING PRESSURE  
Provides a cycle of loading and off-loading by alternating the inflation and deflation of the air bladders.
- Adjustable comfort control settings
- Auto firm
- Power failure alarm and indicator
- CPR hose disconnects for rapid deflation



### OPTIONS

- Available for rent

For more information, please reach out to your Baxter representative.

# SYNERGY AIR ELITE THERAPY SURFACE



## TECHNICAL SPECIFICATIONS

Length (inflated)	79.9" (203 cm)
Width (inflated)	36" (91 cm), 39" (99 cm) or 48" (122 cm)
Height (inflated)	9.4" (24 cm)
Surface Weight	33.1 lb (15 kg)
Air Supply Unit (L x W x H)	16.5" x 6.7" x 10.6" (42 cm x 17 cm x 27 cm)
Air Supply Unit Weight	13.2 lb (6 kg)
Power Requirement	110 V to 120 V, 60 Hz, 5 A
Patient Weight Range	
36" (91 cm) width	86-600 lb (39-272 kg)
39" (99 cm) width	86-1000 lb (39-454 kg)
48" (122 cm) width	100-1000 lb (45-454 kg)

The selection of support surface must be based on a full assessment of the patient's condition and needs, and not solely on the patient's risk for pressure injury.<sup>2</sup>

For safe and proper use of the products mentioned herein, please refer to the appropriate Operator's Manual or Instructions for Use. Baxter reserves the right to make changes without notice in design, specifications and models. The only warranty Baxter makes is the express written warranty extended on the sale or rental of its products.

### References

1. McNichol L, Watts C, Mackey D, Beitz JM, Gray M. Identifying the Right Surface for the Right Patient at the Right Time: Generation and Content Validation of an Algorithm for Support Surface Selection. *Wound Ostomy Continence Nurs.* 2015;42(1):19-37.
2. McNichol L, Mackey D, et al. Choosing a support surface for pressure injury prevention and treatment. *Nursing* 2020; 50(2):41-44

### [bax.to/rental](https://www.baxter.com/rental)

Baxter and Synergy are trademarks of Baxter International Inc. or its subsidiaries.

US-CS158-240002 (v1.0) 12/2024