



Shown with NEW Optional 660lbs Scale

## RESIDENT FLOOR LIFTS

Surprisingly affordable! The new DRIVE Floor Lifts give you more for your money: best-in-class warranty with innovative features and a contemporary, stylish design. Why settle for less when you can have more for about the same cost? Upgrade to Drive DeVilbiss Healthcare today.

### STANDARD FEATURES

- Contemporary soft curve design and colors blend elegantly into any Long Term Care or Assisted Living facility.
- Available in Powered and Manual Base
- Lift control station and battery are smaller and complement design of the lift
- Smart & Safety Technology
  - Redundant lift controls and visual battery fuel indicator can be found on both the Hand Pendant and Lift Control Station
  - Auditable alert battery fuel indicator gently reminds caregiver that the battery is running low
  - O.O.S. (Object Obstruction Sensor) stops lift if resident or caregiver accidentally gets in the way of the boom arm.
  - Audible Alert if the lift detects weight overload
  - Maintenance monitoring:
    - Lift Control Station displays lift utilization and service data. Will also show usage history and has a service reminder to keep you on track for compliance.
    - Utilizing data is key in facilities who look for solutions in meeting compliance with their Safe Handling Programs
  - Stocked available for quick ship

### OPTIONAL FEATURE

Item # FLSCALE  
Resident Floor Lift Scale, Weight Capacity 660 lbs., 300 kg

### SPECIFICATIONS

Boom Height	26" - 76"
Base Width	29" (open)   43" (closed)
Base Height	5.3"
Casters	3" Front, 4.5" Rear
Weight	115 lbs.
Carton Shipping Weight	133 lbs.
Weight Limit	500 lbs   600 lbs.
Warranty	Lifetime on Frame, 5 Years on Electronics

ITEM #	Product	Power/Manual	Wt Cap
FLNP500	Levantar	Manual	500 lbs
FLP500	Levantar	Power	500 lbs
FLNP600	Gravis	Manual	600 lbs
FLP600	Gravis	Power	600 lbs

To Learn More/Request a Demo: Call our LTC Customer Service Hotline at 1-877-374-8763, Contact your Sales Representative, or Visit [DriveMedical.com](http://DriveMedical.com) TODAY!

